

## EXAM GLOVES CATALOGUE

SINGAPORE



# THE WORRY-FREE CHOICE FOR EVERYONE, EVERYWHERE

Your hands are essential to the care you provide to your patients. With 5–15% of healthcare workers susceptible to latex allergy, caring for you and your hands by reducing the exposure to latex reduces the risks and harmful effects of an allergic reaction<sup>1</sup>.

Halyard has created a patented, natural rubber latex-free glove to provide you with the comfort, fit and sensitivity of latex, without the added risk of exposure.

It's the one glove for all your non-surgical needs, facility-wide<sup>2</sup>.



## INNOVATION IN PACKAGING TO REDUCE CROSS CONTAMINATION AND WASTE

As part of our commitment to reducing the risk of HAIs, Halyard introduced the SmartPULL\* glove dispensing system to our STERLING\* and LAVENDER\* Nitrile Exam Glove boxes.

Our advanced packing technology has also reduced the waste in storage space with the larger, 200–250 count per box providing you with more for less.



The smaller tab reduces the risk of contamination and waste by limiting the contact and potential fall-out from adjacent gloves.



The larger tab enables easier access as the box empties.

#### REFERENCES:

- 2010, The Australasian Society of Clinical Immunology and Allergy (ASCIA) http://www.allergy.org.au/health-professionals/hp-information/asthma-and-allergy/latex-allergy (Viewed 26/02/15)
- 2. Halyard offers a line of gloves specifically for use with chemotherapeutic agents.

#### **NITRILE EXAM GLOVES**

#### **STERLING\* NITRILE EXAM GLOVES**

- Powder-Free
- Natural Rubber Latex-Free
- Non-Sterile
- Ambidextrous
- Beaded Cuff
- Single Use Only
- CE Marked
  - Enhanced Wet and Dry Grip • Tested Against 52 Chemotherapy Agents

• Length: 9.5" (241mm)

• Silver









Code	Description	Packaging
50705	STERLING* Nitrile Exam Gloves, X-Small	Eaches/Box = 200 Boxes/Case = 10
50706	STERLING* Nitrile Exam Gloves, Small	Eaches/Box = 200 Boxes/Case = 10
50707	STERLING* Nitrile Exam Gloves, Medium	Eaches/Box = 200 Boxes/Case = 10
50708	STERLING* Nitrile Exam Gloves, Large	Eaches/Box = 200 Boxes/Case = 10
50709	STERLING* Nitrile Exam Gloves, X-Large	Eaches/Box = 170 Boxes/Case = 10

#### **STERLING\* NITRILE-XTRA EXAM GLOVES**

- Powder-Free
- Natural Rubber Latex-Free
- Non-Sterile
- Ambidextrous • Beaded Cuff
- Single Use Only
- Length: 12" (305mm)
- Silver
- Enhanced Wet and Dry Grip
- Tested Against 52 Chemotherapy Agents







Code	Description	Packaging
53137	STERLING* NITRILE-XTRA Exam Gloves, X-Small	Eaches/Box = 100 Boxes/Case = 10
53138	STERLING* NITRILE-XTRA Exam Gloves, Small	Eaches/Box = 100 Boxes/Case = 10
53139	STERLING* NITRILE-XTRA Exam Gloves, Medium	Eaches/Box = 100 Boxes/Case = 10
53140	STERLING* NITRILE-XTRA Exam Gloves, Large	Eaches/Box = 100 Boxes/Case = 10
53141	STERLING* NITRILE-XTRA Exam Gloves, X-Large	Eaches/Box = 100 Boxes/Case = 10

### **PURPLE NITRILE\* EXAM GLOVES**

- Powder-Free
- Natural Rubber Latex-Free
- Non-Sterile
- Ambidextrous • Beaded Cuff
- Single Use Only
- Length: 9.5" (241mm)
- Purple
- CE Marked
- Enhanced Wet and Dry Grip
- Tested Against 29 Chemotherapy









Code	Description	Packaging	
55080	PURPLE NITRILE* Exam Gloves, X-Small	Eaches/Box = 100 Boxes/Case = 10	
55081	PURPLE NITRILE* Exam Gloves, Small	Eaches/Box = 100 Boxes/Case = 10	
55082	PURPLE NITRILE* Exam Gloves, Medium	Eaches/Box = 100 Boxes/Case = 10	
55083	PURPLE NITRILE* Exam Gloves, Large	Eaches/Box = 100 Boxes/Case = 10	
55084	PURPLE NITRILE* Exam Gloves, X-Large	Eaches/Box = 90 Boxes/Case = 10	

#### **PURPLE NITRILE-XTRA\* HIGH RISK EXAM GLOVES**

- Powder-Free
- Natural Rubber Latex-Free
- Non-Sterile
- Ambidextrous
- Beaded Cuff
- Single Use Only
- Purple

- Length: 12" (305mm)Longer Cuff for Extra Protection
- CE Marked
- Enhanced Wet and Dry Grip
- Tested Against 52 Chemotherapy Agents









Code	Description	Packaging
50600	PURPLE NITRILE-XTRA* Exam Gloves, X-Small	Eaches/Box = 50 Boxes/Case = 10
50601	PURPLE NITRILE-XTRA* Exam Gloves, Small	Eaches/Box = 50 Boxes/Case = 10
50602	PURPLE NITRILE-XTRA* Exam Gloves, Medium	Eaches/Box = 50 Boxes/Case = 10
50603	PURPLE NITRILE-XTRA* Exam Gloves, Large	Eaches/Box = 50 Boxes/Case = 10
50604	PURPLE NITRILE-XTRA* Exam Gloves, X-Large	Eaches/Box = 50 Boxes/Case = 10

## PURPLE NITRILE\* MAX EXAM GLOVES

- Powder-Free
- Natural Rubber Latex-Free
- Non-Sterile
- Ambidextrous
- Beaded Cuff
- Single Use Only
- Purple

- Longer Length: 16 inches¹ (406mm)
- CE Marked
- Enhanced Wet and Dry Grip
- Heavier weight: 20 grams<sup>2</sup>
- Increased fingertip thickness
- Tested Against 52 Chemotherapy Agents and 25 Commonly Used Chemical and Detergents









Code	Description	Packaging
44992	PURPLE NITRILE* MAX Exam Gloves, Small	Eaches/Box = 50 Boxes/Case = 8
44993	PURPLE NITRILE* MAX Exam Gloves, Medium	Eaches/Box = 50 Boxes/Case = 8
44994	PURPLE NITRILE* MAX Exam Gloves, Large	Eaches/Box = 50 Boxes/Case = 8
44995	PURPLE NITRILE* MAX Exam Gloves, X-Large	Eaches/Box = 50 Boxes/Case = 8

## LAVENDER\* NITRILE EXAM GLOVES

- Powder-Free
- Natural Rubber Latex-Free
- Non-Sterile
- Ambidextrous
- Beaded Cuff
- Single Use Only
- Length: 9.5" (241mm)
- Lavender
- Enhanced Wet and Dry Grip
- Tested Against 43 Chemotherapy Agents





Code	Description	Packaging
52816	LAVENDER* Nitrile Exam Gloves, X-Small	Eaches/Box = 250 Boxes/Case = 10
52817	LAVENDER* Nitrile Exam Gloves, Small	Eaches/Box = 250 Boxes/Case = 10
52818	LAVENDER* Nitrile Exam Gloves, Medium	Eaches/Box = 250 Boxes/Case = 10
52819	LAVENDER* Nitrile Exam Gloves, Large	Eaches/Box = 250 Boxes/Case = 10
52820	LAVENDER* Nitrile Exam Gloves, X-Large	Eaches/Box = 230 Boxes/Case = 10

#### REFERENCES:

- 1. Specification target. Data on file
- X-Large size

## AQUASOFT\* NITRILE EXAM GLOVES

- Powder-Free
- Natural Rubber Latex-Free
- Non-Sterile
- Ambidextrous
- Beaded Cuff
- Single Use Only
- Length: 9.5" (241mm)
- CE Marked
- Textured fingertips
- Tested Against 34 Chemotherapy Agents
- Higher count per box saves space and waste: 3000 units per case / 300 gloves per inner









Code	Description	Packaging
43932	AQUASOFT* Nitrile Exam Gloves, X-Small	Eaches/Box = 300 Boxes/Case = 10
43933	AQUASOFT* Nitrile Exam Gloves, Small	Eaches/Box = 300 Boxes/Case = 10
43934	AQUASOFT* Nitrile Exam Gloves, Medium	Eaches/Box = 300 Boxes/Case = 10
43935	AQUASOFT* Nitrile Exam Gloves, Large	Eaches/Box = 300 Boxes/Case = 10
43936	AQUASOFT* Nitrile Exam Gloves, X-Large	Eaches/Box = 250 Boxes/Case = 10

## **GLOVE DISPENSERS**



#### **HALYARD\* EVERYDAY GLOVE DISPENSERS**

Code	Description	Packaging
50201	Exam Glove Dispenser (Suits 9.5" - 12" gloves), Clear Plastic, Single, 96mm (D) x 144mm (W) x 253mm (H)	Eaches = 1
50202	Exam Glove Dispenser (Suits 9.5" - 12" gloves), Clear Plastic, Double, 96mm (D) x 274mm (W) x 256mm (H)	Eaches = 1
50203	Exam Glove Dispenser (Suits 9.5" - 12" gloves), Clear Plastic, Triple, 96mm (D) x 406mm (W) x 256mm (H)	Eaches = 1

### **HALYARD\* PURPLE NITRILE\* MAX GLOVE DISPENSERS**





Code	Description	Packaging
47467	Exam Glove Dispenser Extra Long (Suits 16" gloves), Clear Plastic, Single, 75mm (D) x 139mm (W) x 409mm (H)	Eaches = 1
47466	Exam Glove Dispenser Extra Long (Suits 16" gloves), Clear Plastic, Triple, 69 mm (D) x 399mm (W) x 409 mm (H)	Eaches = 1



Infection Prevention Product Solutions Halyard\* Education Foundation Tools & Best Practice Commitment to Excellence and Sustainability **Experienced Customer Support** Knowledgeable Sales Force Innovation & Technology



#### O&M Halyard Singapore Pte Ltd

United Square, 101 Thomson Road #25-05 S(307591) Tel: +65 8030 5277

www.halyardhealth.sg









Halyard makes no representation that its products described in this catalogue will provide complete protection from exposure to contaminants, fluids or micro-organisms or will eliminate the risk of contracting infectious diseases.